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English

Device Description

Silverlon® Flex Wound Contact Negative Pressure Dressings are sterile, single layer, non-adherent, antimicrobial barrier dressings, intended for topical use in wound management and in negative pressure wound therapy (NPWT), to reduce risk of wound infection.

Available as:

Code	Size	Code	Size
WCD22-NPS(-FLX)	5 x 5 cm	WCD48-NPS(-FLX)	10 x 20 cm
WCD44-NPS(-FLX)	10 x 10 cm	WCD88-NPS(-FLX)	20 x 20 cm

Silverlon® Flex Wound Contact Negative Pressure Dressings contain silver to provide an ancillary, local antimicrobial effect to reduce the risk of wound infection by organisms sensitive to silver. Ionic silver is a broad spectrum antimicrobial, which has been shown to be effective against a range of pathogens associated with wound infection. Silverlon® dressings have been tested in vitro and found effective against microorganisms such as: Enterococcus faecalis, Staphylococcus aureus, Pseudomonas aeruginosa, Acinetobacter Baumannii, Enterobacter cloacae, Staphylococcus epidermidis and Klebsiella pneumoniae. Silverlon® Flex Wound Contact Negative Pressure Dressings may be used with all types of negative pressure wound therapy porous foam/gauze dressings and equipment.

In vitro testing has shown Silverlon® Flex Wound Contact Negative Pressure Dressings to be MRI-compatible.

Silverlon® Flex Wound Contact Negative Pressure Dressings are recommended for use by professional healthcare providers and clinical practitioners.

Intended Purpose

Silverlon® Flex Wound Contact Negative Pressure Dressings are sterile, single layer, non-adherent, antimicrobial barrier dressings, intended for topical use in wound management and in negative pressure wound therapy (NPWT), to reduce risk of wound infection.

Indications

Silverlon® Flex Wound Contact Negative Pressure Dressings are professional use wound dressings indicated for use on acute and chronic wounds and in conjunction with negative pressure wound therapy (NPWT), for local management of:

- 1st and 2nd degree burns
- Donor and graft site wounds
- Diabetic, pressure and vascular ulcers
- Traumatic and surgical wounds
- Infected wounds*

* Silverlon® Flex Wound Contact Negative Pressure Dressings may be used on those infected wounds listed above. Where used on infected wounds, the infection should be treated as per local clinical protocol.

Silverlon® Flex Wound Contact Negative Pressure Dressings are indicated for a total duration of contact (i.e. where treatment involves consecutive application of individual dressings) of up to 30 days.

Warnings

- Do not use past expiration date on the product packaging.
- Do not use if pouch is damaged or open.
- Do not use petroleum-based ointments or creams under Silverlon® Flex Wound Contact Negative Pressure Dressings.
- Avoid contact with electrodes or conductive gels during electronic measurements e.g. EEG and ECG.
- Do not moisten Silverlon® Flex Wound Contact Negative Pressure Dressings with hydrogen peroxide or povidone iodine.
- Silverlon® Flex Wound Contact Negative Pressure Dressings dressings should not be used for children, babies, premature infants, neonates, pregnant and lactating women due to limited clinical data.
- If patients experience any allergic reactions such as blistering, itching, irritation, redness, skin rashes or discoloration, discontinue use and/or contact healthcare providers.
- In the event of clinical infection, topical silver does not replace the need for systemic therapy or other adequate infection treatment.
- Silverlon® Flex Wound Contact Negative Pressure Dressings are intended for single use only and should not be re-used. Re-use may adversely affect the dressing's performance characteristics and also presents a risk of infection to the patient.
- Total duration of use of Silverlon® Flex Wound Contact Negative Pressure Dressings (i.e. where treatment involves consecutive application of individual dressings) should not exceed 30 days.
- Silverlon® Flex Wound Contact Negative Pressure Dressings are not intended for direct application to exposed underlying structures (i.e. blood vessels, bone, tendons, nerves).

Silverlon®

Flex Wound Contact Negative Pressure Dressings

Publication Date: 06/2025

Emurgo Europe
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands
Tel. +(31) (0) 70 345-8570

Argentum Medical LLC dba Bravida Medical
2571 Kaneville Ct
Geneva, IL 60134 USA
Tel. (+1) 888.551.0188
Fax. (+1) 888.558.9923
www.bravidamedical.com

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