Emergo Europe
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands
Tel. +(31) (0) 70 345-8570

Argentum Medical, LLC dba Bravida Medical
2571 Kaneville Ct
Geneva, IL 60134 USA
Tel. (+1) 888.551.0188
Fax. (+1) 888.558.9923
www.bravidamedical.com

English

FreeDerm® Adhesive Remover Instructions for Use

INGREDIENTS: Siloxanes

WARNING! Highly Flammable
Avoid contact with eyes and mouth.
For external use only. Avoid inhalation.

NOT MADE WITH NATURAL RUBBER LATEX

- Removes hydrocolloid on contact
- Helps remove tapes and bandages
- Quick drying – no residue

If redness or irritation occurs, discontinue use and consult your healthcare professional.

DO NOT FLUSH

Any serious incident or complaint in relation to the device report to Manufacturer Argentum Medical, LLC dba Bravida Medical, 2571 Kaneville Court, Geneva, IL 60134 Phone no. + 1 888.551.0188 and To EC Rep Emergo Europe Westervoortsedijk 60; 6827 AT Arnhem The Netherlands Tel. +(31) (0) 70 345-8570 and report to local competent authority.

©Argentum Medical, LLC dba Bravida Medical. All rights reserved. FreeDerm® is a registered trademark of Argentum Medical, LLC dba Bravida Medical. Made in the USA.
Patents: <https://bioderminc.com/patents/>.

HOW TO USE: Instantly removes hydrocolloid, dressing tapes and other adhesives without irritation. Flash dries in 15 seconds without leaving a residue which allows for immediate reapplication.

SPRAY BOTTLE DIRECTIONS: For external use only on normal, intact skin. Spray edges of adhesive until wet. Gently lift the adhesive edge from the skin. While lifting, spray the exposed adhesive until completely released.

WIPE DIRECTIONS: For external use only on normal, intact skin.

1. Hold foil packet vertically, tear at notch to open and remove cloth from packet. Keep cloth folded to retain maximum saturation. Keep the packet held upright to avoid spilling the excess FreeDerm solution.
2. Drip the excess FreeDerm from the packet onto the top edge of the dressing, allowing the solution to drip downwards to start releasing the adhesive bond between the skin and dressing.
3. Once the dressing edge is lifted, use the folded cloth to wipe underneath the dressing as you gently pull the dressing off.
4. Once the dressing is completely removed, wipe the area to remove any residue adhesive.

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