

CathGrip®

Double Strap

LC-BOM51300BX-01 A



Manufactured by: Argentum Medical, LLC 2571 Kaneville Ct.,
Geneva, IL 60134 | www.BravidaMedical.com
Made in USA with global components.

bravida
medical

©Argentum Medical, LLC DBA Bravida Medical. All rights reserved.
BioDerm® and CathGrip® are registered trademarks of BioDerm,
Inc., DBA Bravida Medical. BioPlus+™ is a trademark of BioDerm,
Inc., DBA Bravida Medical. Patents: <https://bioderminc.com/patents/>.



STERILE R

10°C
(50°F)



30°C
(85°F)

Not Made with Natural Rubber Latex

INSTRUCTIONS FOR USE - A,B,C,D

- A - Apply Alcohol Wipe on Application Site. Allow to Dry 15 Seconds
- B - Apply BioPlus+ Wipe on Application Site. Allow to Dry 15 Seconds
- C - Apply CathGrip on Application Site. Press CathGrip on the Skin for 15 Seconds
- D - Place Tube inside Strap System. Pull and Secure Strap(s)

Pull Open Here